




MAX WOCHER & SON, CINCINNATI, O.



D/..1050. Fig. B



D/..1050. Fig. A



D/..1061

D/.. 978. **Dommer's** urethral and rectal electrodes for faradization of the posterior urethra, prostate, etc. per set \$ 8 00


D/..1050. **Bierhoff-Freudenberg's** cystoscopic prostate incisor, with cable cords and extra lamp, in case. 125 00

Fig a. Showing instrument used as incisor

" b. " " " cystoscope.

This instrument enables the operator to examine before and after the cut, and as the cystoscope can be readily withdrawn, without removing the entire instrument, the bladder can be irrigated, refilled, etc., whenever desirable and the cystoscope re-inserted through the tube.

D/..1061. **Transformer** for 52 or 104 volt alternating current, specially wound for the Rottini operation. 30 00



D/..978

ALL INSTRUMENTS BEARING THE K.-S. CO.'S TRADE MARK ARE SOLD UNDER GUARANTEE



MALE URETHRAL INSTRUMENTS.