


ELECTRICAL BATTERIES.



THE NEW No. 4 DRY CELL MEDICAL APPARATUS.

No ACID. No LIQUIDS.
No REPAIRS.

Always Ready, Neat, Clean,
Compact, Portable.

This battery is divided into two sections. The upper receptacle holds the following: Two Metal Tube Electrodes, two Sponge Electrodes, two Insulating Handles, for use with either Hand Electrodes or Sponge Electrodes, one Foot Plate and one pair Conducting Cords. The lower receptacle, as shown on the illustration, contains the Battery and Coil.

Price to Patients, \$7 50

To Physicians, 4 50

Price for extra Cells, each, 25

Three different currents are obtained from the apparatus: The mildest or primary current (alone) by placing the cord tips in the contact holes marked 1 and 2. By placing the tips in 2 and 3, the secondary current (alone) is obtained, and by placing them in 1 and 3, the combined currents are given out.

The strength of all the currents can be varied by moving the regulator or core-shield of the coil; the further out this shield is drawn the stronger the currents become.

The Dry Battery which is used with this apparatus, being a standard staple commercial article, made only in one size and shape, is easily obtainable at any Electrical Supply store, or wherever this No. 4 apparatus is on sale.

One cell will last from four to twelve months, according to the length and frequency of applications.

DRY CELLS

For Faradic or Galvanic Batteries, Cabinets and Electric Bells.

Price.....25 cts. each.

Special discount in quantities.

All Instruments designated with a * are illustrated.

