


MAX WOCHER & SON, CINCINNATI, O.



ELECTROLYTIC INSTRUMENTS FOR THE URETHRA.

D/..868.	Casper's electrolytic needle, for the treatment of hypertrophia prostatae	4 50	D/..898.	Lang's electrolytic stricture dilating sound, straight and short, for the anterior parts.:	\$ 6 00
D/..874.	Kollmann's electrolytic needle, blunt pointed	3 00	D/..904.	" flexible stricture sound electrode. . .	9 00
D/..880.	" " " sharp "	3 00	D/..910.	Newman's electrolytic sound, straight, hard rubber cover	1 00
The points of electrolytic needles are of platinum. The negative pole is attached to handle while the positive pole is attached to an electrode applied to the thigh; strength of the current required, 4 to 6 milliamperes.			D/..916.	" " " curved, hard rubber cover	1 50
D/..886.	Oberlaender's fork needle.....	6 00	D/..922.	" " " curved, canulated, with filiform bougie	5 00
With this instrument either both negative poles or one positive and one negative pole can be applied. In the latter case only a weak current is needed.			D/..928.	Fort's electrolytic bougie, with filiform guide for linear electrolysis.	5 00
D/..892.	Lang's electrolytic stricture dilating sound, curved....	6 00	D/..934.	" " " with filiform guide to unscrew.....	6 00

ALL INSTRUMENTS BEARING THE K.-S. CO.'S



TRADE MARK ARE SOLD UNDER GUARANTEE.

MALE URETHRAL INSTRUMENTS.