


# MAX WOCHER & SON, CINCINNATI, O.



## ELECTROLYTIC INSTRUMENTS FOR THE URETHRA.—Continued.

- D/..940. Kollmann's electrode for galvanization or electrolysis of the mucous membranes of the urethra..... \$5 00  
 The instrument consists of a hard rubber covered silver tube with perforations on the side running all through. When introduced, the olive tip on end prevents the fluid from passing into the bladder. The conical part on other end is to close the orifice and prevent exit of fluid. When urethra is filled electricity is applied.

## URETHRAL ELECTRODES.

- |           |                                                                                             |      |
|-----------|---------------------------------------------------------------------------------------------|------|
| D/..952.  | Kollmann's intra urethral electrode, for application through endoscope.....                 | 1 50 |
| D/..958.  | Intra urethral electrode, can be applied through electro-urethroscope when illuminated..... | 1 50 |
| D/..964.  | Bladder electrode, flexible, with insulated covering.....                                   | 1 75 |
| D/..970.  | " " with water canal for filling the bladder.....                                           | 5 00 |
| D/..982.  | " " bipolar.....                                                                            | 3 00 |
| D/..988.  | " " for female urethra.....                                                                 | 1 75 |
| D/..994.  | Stembo's perineal electrode.....                                                            | 4 00 |
| D/..1000. | Zuelzer's scrotal ".....                                                                    | 2 50 |

ALL INSTRUMENTS BEARING THE K.-S. CO.'S  TRADE MARK ARE SOLD UNDER GUARANTEE.

MALE URETHRAL INSTRUMENTS.